

# Pavlik Harness

*A soft bracing system for the  
treatment of Hip Dysplasia*



## The DM Orthotic Pavlik Harness:-

- A specially designed soft brace
- Keeps hips flexed and abducted
- Keeps femoral head in the correct position
- Encourages bones to grow normally

Hip dysplasia is used to describe a range of conditions from a hip that is mildly dysplastic, centrally located and stable, to one that is severely dysplastic and dislocated.

In children younger than 6 months the conventional treatment of dysplasia is bracing.

### ***The DM Orthotic Pavlik Harness can be worn:-***

- from birth where the hip is unstable
- in older infants where there is subluxation or actual dislocation of the hip

The degree of flexion and abduction can be adjusted according to the specific needs of the child.

The colour coded straps help for ease of re-application.

The DM Orthotic Pavlik Harness is available with or without booties.

### ***The DM Orthotic Pavlik Harness is available in the following sizes:-***

	Chest Circumference
X-Small	25 - 34 cm
Small	34 - 43 cm
Medium	43 - 52 cm
Large	51 - 61 cm

**For more information please contact +44 (0) 1209 219205  
[www.dmorthotics.com](http://www.dmorthotics.com) [admin@dmorthotics.com](mailto:admin@dmorthotics.com)**